



Tvinnargatan 25 · SE-507 30 Brämhult Tfn: 033-22 56 30 · info@elastocon.se

www.elastocon.se

"The device is working perfectly."

Dr. Christian Heyde
Intelligence Engine
adidas



Quantify the warm feel and cool touch of textiles and fabrics via ASTM D7984

## CHARACTERIZING COOL TOUCH


## via ASTM D7984 with C-Therm Tx

Touch is one of our critical senses in perceiving the world – everything from the clothes we wear to the bedding we sleep in. Human skin is very good at detecting differences in a material's ability to transfer heat, such as the warmth of a fleece sweater compared to the coolness of leather. This material property is known as thermal effusivity, and as a metric, it can be used to quantify a textile's ability to exchange thermal energy between skin and fabric.

## Tx Touch Experience Platform

C-Therm's Tx Platform allows for rapid, repeatable and controlled measurement of thermal effusivity and thermal conductivity of materials. Employing the Modified Transient Plane Source Sensor with a compression test apparatus, a limited (~1 °C) heat pulse is emitted over a short test time (~1s) and the voltage drop on the sensor surface is precisely monitored. This allows the thermal properties of the material in contact with the sensor to be determined.





## The Warm Feel - Cool Touch Product Performance Index



Modified Transient Plane Source Sensor

Units of Thermal Effusivity are W s1/2 / m2K